

**Departure Instructions**  
**UNOLS East Coast Van Pool**  
**Dry Lab #2**  
**Sonic Barriers #625.2.01-1**

The following list is intended to aid in removing the lab from the ship and making it ready to return to the UECVP. Any questions, concerns, or discrepancies noted should be referred to the UECVP Manager (Timothy W. Deering 302-645-4338 or [deering@udel.edu](mailto:deering@udel.edu)).

- The lab should be clean and dry including the inside of the hood, freezer and refrigerator. No scientific equipment, chemicals, supplies, or waste are to be left in the lab.
- Complete and return the Evaluation Form in this packet.
- Turn off all equipment including the refrigerator, freezer, and HVAC unit(s). Unplug the water heater.
- Disconnect the fresh water, salt water, and sink drain hose. Drain and stow them in the lab. Replace the cam-lock covers and caps on all fittings.
- Remove the hood gooseneck and any cable pass cowls and stow in the lab on the racks provided. Put the cover plates on all wire passes.
- Put cover plates on all windows.
- Turn the switch on the emergency light to the “Off” position.
- Disconnect the electrical cable and stow inside the lab. Place the cover on the electrical receptacle.
- Make sure all drawers are latched shut, loose items are secured, and the refrigerator and freezer doors are latched shut.
- Lock the lab and make suitable arrangements to return the keys to the UECVP (send them with the truck driver, mail them, etc.).