

UNOLS Isotope Van LSC Start-Up and Shut-Down Procedure

Start-Up

1. Start the LSC first—ON/OFF switch is located on the lower right hand side towards the front of the unit.
2. Start-up the lap top, there is not a password.
3. Start-up the printer--- ON/OFF switch is located on the back left hand side.

Shut Down

1. Turn off printer
2. Turn off laptop
3. Turn off LSC; remember to make sure everything is removed from inside the LSC before locking up the van for shipping.

If the LSC starts up when you power up the van shut everything down. Then follow the Start-Up Procedure.